

# Capp *Bravo*

## SINGLE CHANNEL PIPETTES

Manual pipetting has never been easier! **Capp Denmark** has combined lightweight ergonomic design with the latest innovations in manual liquid handling to create **Capp *Bravo***, the most user friendly manual pipette on the market today.

**Capp *Bravo*** is constructed of the highest quality autoclavable materials and features feather light plunger and ejection forces. There is a built-in oversized indicator for easy volume monitoring and precise adjustment. With **Capp *Bravo*** all of these features are packaged in a modern, streamlined, world-renowned Danish design.

Experience Capp's symphony of precise measurements and see for yourself why it is called Bravo!

### Safety:

- Fully Autoclavable
- Positive indicator lock

### Ergonomics:

- Ergonomically constructed hook for comfortable handling
- Smooth tip mounting and ejection

### Controls:

- Oversized indicator
- Reduced force feather touch plunger and ejector

### Unique Features:

- Metal shaft with double o-rings for ultra micro sizes
- Ultra lightweight construction



# Ordering Info

## Capp *Bravo* Starter Packs

| Cat. No. | Description  |
|----------|--|
| SK-03-BB | Three pipettes free of choice, C-05 stand, 3 racks of tips and 1 timer |

## Pipette Stands

| Cat. No. | Description                        |
|----------|------------------------------------|
| C-05     | Pipette stand for up to 5 pipettes |



## Capp *Bravo* variable volume

| Cat. No. | Description   | Inaccuracy % | Imprecision % | Colour |
|----------|---|--------------|---------------|--------|
| B02-1    | Capp pipette, variable vol. 0.1-2 $\mu$ l, Bravo    | 12.00 / 2.00 | 6.00 / 1.00   | ●      |
| B10-1    | Capp pipette, variable vol. 0.5-10 $\mu$ l, Bravo   | 2.50 / 1.00  | 1.80 / 0.50   | ○      |
| B20-1    | Capp pipette, variable vol. 2-20 $\mu$ l, Bravo     | 5.00 / 1.00  | 1.50 / 0.30   | ●      |
| B50-1    | Capp pipette, variable vol. 5-50 $\mu$ l, Bravo     | 2.00 / 1.00  | 1.20 / 0.30   | ●      |
| B100-1   | Capp pipette, variable vol. 10-100 $\mu$ l, Bravo   | 2.50 / 0.70  | 1.00 / 0.20   | ●      |
| B200-1   | Capp pipette, variable vol. 20-200 $\mu$ l, Bravo   | 2.50 / 0.60  | 0.70 / 0.20   | ●      |
| B1000-1  | Capp pipette, variable vol. 100-1000 $\mu$ l, Bravo | 2.00 / 0.80  | 0.60 / 0.20   | ●      |
| B5000-1  | Capp pipette, variable vol. 1-5mL, Bravo            | 1.20 / 0.60  | 0.30 / 0.15   | ●      |
| B10000-1 | Capp pipette, variable vol. 1-10mL, Bravo           | 3.00 / 0.60  | 0.60 / 0.15   | ●      |

## Capp *Bravo* fixed volume *(customized volumes available on request)*

| Cat. No.  | Description                                  | Inaccuracy % | Imprecision % |
|-----------|--|--------------|---------------|
| B01-1F    | Capp pipette, fixed vol. 1 $\mu$ l, Bravo    | 2.50         | 1.80          |
| B02-1F    | Capp pipette, fixed vol. 2 $\mu$ l, Bravo    | 2.00         | 1.00          |
| B05-1F    | Capp pipette, fixed vol. 5 $\mu$ l, Bravo    | 1.20         | 0.80          |
| B10-1F    | Capp pipette, fixed vol. 10 $\mu$ l, Bravo   | 1.00         | 0.50          |
| B20-1F    | Capp pipette, fixed vol. 20 $\mu$ l, Bravo   | 1.00         | 0.30          |
| B50-1F    | Capp pipette, fixed vol. 50 $\mu$ l, Bravo   | 1.00         | 0.30          |
| B100-1F   | Capp pipette, fixed vol. 100 $\mu$ l, Bravo  | 0.70         | 0.20          |
| B200-1F   | Capp pipette, fixed vol. 200 $\mu$ l, Bravo  | 0.60         | 0.20          |
| B500-1F   | Capp pipette, fixed vol. 500 $\mu$ l, Bravo  | 0.80         | 0.15          |
| B1000-1F  | Capp pipette, fixed vol. 1000 $\mu$ l, Bravo | 0.80         | 0.20          |
| B5000-1F  | Capp pipette, fixed vol. 5mL, Bravo          | 0.60         | 0.15          |
| B10000-1F | Capp pipette, fixed vol. 10mL, Bravo         | 0.60         | 0.15          |