

Refrigerated Type (Medium Capacity)



CN-350R



CN-350RDS



CN-350RD

特點：

- 採用無刷馬達。
- LCD數位面板，具溫度、轉速及離心力雙重切換、設定及顯示功能。
- 10組記憶系統，可預先設定常用參數。
- 不平衡斷電、過載、故障指示安全防護裝置。
- 具預冷功能。

Feature :

- Brushless Motor
- LCD digital panel with temperature, RPM & RCF switch, setting and display functions
- 10 memory programs for pre-setting parameters
- Auto power off for rotor imbalance, overload, error indication and safety protection device
- With the pre-cooling Function



掀蓋內部 Inside

| Model | CN-350R | CN-350RD | CN-350RDS |
|-------------------------|----------------------|-----------------|-----------------|
| 轉速 Speed (rpm) | 500~6,000 | | |
| 最大離心力 Max. RCF (xg) | 4,270 | | |
| 最大容量 Max. Capacity (ml) | 400 | | |
| 溫度範圍 Temperature range | 0°C~40°C | | |
| 定時器 Timer (min:sec) | 99 : 59 | | |
| 重量 Weight (kg) | 85 | 55 | 60 |
| 尺寸 Size (W x L x H mm) | 430 x 580 x 850 | 650 x 550 x 350 | 440 x 720 x 350 |
| 電源 Power | AC 110/220V, 50/60Hz | | |

● 離心試管另購 · Tubes not included.

■ Rotor sets for CN-350R, CN-350RD, CN-350RDS



| Rotor | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|-------|---------------------|----------------------|------------------|---------------|
| RS-50 | 1.5ml x 48 | 9.5~10.5 x 38~50 | 3,500 | 1,320 |
| | 6ml x 36 | 11~12 x 85~106 | 3,500 | 1,940 |
| | 6ml x 40 | 11~12.5 x 63~87 | 3,500 | 1,630 |
| | 10ml x 24 | 14~15.5 x 87~108 | 3,500 | 1,980 |
| | 15ml x 8 | 15~17.3 x 90~115 | 3,500 | 1,980 |
| | 15ml x 12 (plastic) | 15~17 x 89~115 | 3,500 | 2,000 |
| | 15ml x 16 | 15~17.3 x 90~115 | 3,500 | 1,980 |
| | 50ml x 4 (plastic) | 26.5~29.5 x 90~120 | 4,000 | 2,560 |
| | 50ml x 4 (glass) | 32~35 x 90~112 | 4,000 | 2,490 |
| | 100ml x 4 | 42~45 x 97~112 | 3,500 | 1,980 |



| Rotor | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|-------|---------------------|----------------------|------------------|---------------|
| RS-55 | 1.5ml x 48 | 9.5~10.5 x 38~50 | 3,500 | 1,320 |
| | 6ml x 36 | 11~12 x 85~106 | 3,500 | 1,940 |
| | 6ml x 40 | 11~12.5 x 63~87 | 3,500 | 1,630 |
| | 10ml x 24 | 14~15.5 x 87~108 | 3,500 | 1,980 |
| | 15ml x 8 | 15~17.3 x 90~120 | 4,500 | 3,270 |
| | 15ml x 12 (plastic) | 15~17 x 89~120 | 3,500 | 2,000 |
| | 15ml x 16 | 15~17.3 x 90~120 | 4,500 | 3,270 |
| | 50ml x 4 (plastic) | 26.5~29.5 x 90~120 | 5,000 | 4,000 |
| | 50ml x 4 (glass) | 32~35 x 90~112 | 5,000 | 3,890 |
| | 100ml x 4 | 42~45 x 97~112 | 4,500 | 3,270 |

■ Rotor sets for CN-350R, CN-350RD, CN-350RDS



| Rotor | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|------------|---------------------|----------------------|------------------|---------------|
| RA-1024 | 10ml x 24 | 15~17.3 x 90~107 | 3,500 | 1,760 |
| RA-1508-1 | 15ml x 8 | 15~17.3 x 92~130 | 6,000 | 3,750 |
| RA-1512H-1 | 15ml x 12 | 15~17.3 x 92~122 | 6,000 | 4,270 |
| RA-1516 | 15ml x 16 | 15~17.3 x 90~115 | 4,000 | 2,220 |



| Rotor | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|---------|---------------------|----------------------|------------------|---------------|
| RA-5004 | 50ml x 4 | 27.5~29.5 x 84~135 | 6,000 | 3,870 |
| RA-5006 | 50ml x 6 | 27.5~29.5 x 84~135 | 6,000 | 3,870 |
| RA-5008 | 50ml x 8 | 27.5~29.5 x 84~116 | 6,000 | 4,270 |