



NEW Water Bath






Features

- Equipped overheat protection sensor.
- Time set function and constant temperature.
- Uniquely designed smooth exterior and mono piece bath.
- Equipped internal bottom for the safety of the user and protection of the unit.

Specifications

Made in Korea

Spec. / Model		WD-11N	WD-23N
Type		Digital	
Controller		PID Controller	
Range		Ambient 5~99 °C	
Temp	Uniformity	±0.5°C	
Size	Internal (W x D x H, mm)	300x240x200	515x330x200
	External (W x D x H, mm)	360x370x280	580x410x280
Safety		Over Temp. Cut-Off	
Material	Internal	Seamless Stainless steel	
	External	Steel, Powder coating	
Electric Requirements		220VAC 50/60Hz	
Heater		1000W	2400W



Model : WD-23N

(mm)

	Model	W	D	H	W1	D1	H1	g
Dimension	WD-11N	360	375	285	303	240	200	33
	WD-23N	580	415	281	498	295	200	33

Internal size