

FW100

High-Speed Universal Disintegrator applies to industries, agriculture, mining enterprises, health and medical units, coal, geological and other scientific research institutions for grain, soil, medicine, dinas and mineral substance grinding treatment.

Characteristics

1. Stainless steel smashing chamber ensures grinding accuracy.
2. Chamber roof cover with hyaline window for convenient observation of grinding results.
3. Characterized by small volume, high efficiency, easy operation, attractive appearance and extensive applications.



Specification

Model	Voltage (V)	Motor Power (W)	Speed (r/min)	Smashing (mesh)	Product Size (mm)	Packing Size (mm)	N. W. (Kg)	G. W. (Kg)
FW100	220	460	24000	60-200	Φ 140× 290	325× 180× 190	3.3± 0.3	3.4

