

Grinder

RT Series

Raw material will be pulverized in short time with extreme high speed RPM (25,000 ~ 30,000 RPM) crushing cutter spinning. RT series is the most suitable for usage of research institutions, laboratory and individual working.



Feature

- The best choice and most popular products which meet the average user demand.
- The machine is suitable and will be the best choice for university, pharmaceutical factory, agriculture and laboratory.
- Crushing chamber made of stainless steel, the installation of safety for safer operation.
- Crushing chamber made of stainless steel, Motor is made of aluminium alloy, light weight, fast heat dissipation, Installation of overload safety switch to prevent coil burnt.

Specification

Model	Electric Power	Voltage	Volume	R.P.M	Crushing Chamber	Dimension (L*W*H mm)	Weight (KG)
RT-02	820W	110~240V	150g	30000	78*113	150*150*300	3.5
RT-04	1200W	110~240V	250g	24000	85*139	180*160*340	6.5
RT-08	1800W	110~240V	350g	24000	140*174	200*170*390	9.6

