



Incubator Shaker

Features


- Accurate temperature controls by Digital PID controller.
- Accurate speed control and digital mark function made for user convenience.
- Circulation device in internal chamber have a little noise and keep overall temperature evenly.
- Overheat, over current prevention device for users and safety function to stop shaking device automatically when door open.
- Various holder and plate available according to users spec.
- Little noise designed structurally

Specifications


Made in Korea

Spec. / Model		SI-64	SI-150
Temp	Range	Amb +5 °C ~ 60 °C	
	Controller	Digital PID Controller	
	Accuracy	±0.2 °C	
Size	Internal (W x D x H, mm)	400x400x400	560x560x500
	External (W x D x H, mm)	470x580x630	630x740x750
Safety	Electrical	Over temp limiter, Door safety	
Material	Internal	Stainless steel(STS #304)	
	External	Steel, Powder coating	
Speed Range		20 ~300 rpm	
Shaking Stroke		20 mm Orbital Motion	
Electric Requirements		220VAC. 50/60Hz	
power consumption		330W	430W

 Distinction



 Structure



Shaking Incubator(SI-100 / SI-100R)



Features

- Accurate temperature controls by Digital PID controller.
- Accurate speed control and digital mark function made for user convenience.
- Circulation device in internal chamber have a little noise and keep overall temperature evenly.
- Overheat, over current prevention device for users and safety function to stop shaking device automatically when door open.
- Various holder and plate available according to users spec.
- Little noise designed structurally.


Specifications

Made in Korea

Spec. / Model		SI-100	SI-100R
Temp	Range	Amb +5℃~60℃	0℃~60℃
	Controller	Digital PID Controller	
	Accuracy	±0.2℃ at 37℃	
Size	Internal (WxDxH, mm)	530x560x380 (112L)	
	External (WxDxH, mm)	590x850x770	590x850x880
Safety	Electrical	Over temp limiter, ELB(Earth Leakage Circuit Breaker), Pause stop	
Material	Internal	Stainless steel (#304)	
	External	Steel, Powder coating	
Speed Range		20 ~300 rpm	
Shaking Stroke		20mm Orbital Motion	
Light System		LED Lamp 4W	
Refrigerator		None	1/4 HP
Electric Requirements		220VAC. 50 or 60Hz	

Shaking Incubator

Distinction



shock absorber

lamp

Spring Wire Rack

PID Controller

* Accurate temp. control and speed control by the touch PID Controller

* Various flasks are available.

* Easy observation during operation.

*Smooth door operation

Structure

