

Digital Friability Tester Model: FT-20



Dual Drums



Front Load



10 Degrees Tilt



25 rpm

Features :

- Complies to USP and other pharmacopoeias
- Compact, Simple & Easy to operate
- Cost Effective
- Supports friability drum with single scoop (Standard), friability drum with double scoop and abrasion drum
- Auto discharge of tablets after friability test with single scoop & Double Scoop friability drum
- Dual mode of operation timer mode or revolution count mode
- Menu driven operation makes user friendly
- Selectable types of operation with Friability Test, Abrasion (Attrition) Test with single or double drum
- 10 deg. tilt facility
- Easy to validate
- Calculation of percentage weight loss and displays the results of selected test type.
- Power failure prompt to discard or continue the test
- Password protection for Real Time Clock**
- IQ/OQ Documents

** Due to continuous product development the image might vary*

Specifications :

- Display : 20 x 4 LCD Display with backlit
- Operation : Through front operating panel with alpha numeric entries
- RPM : 25 Fixed RPM
- Mode : Dual Mode Timer Mode & Count Mode (1 min to 99hrs: 59min: 59sec)
- Accessories : Standard
- Power : 230V AC ~50Hz 1Amp
: 110V AC ~ 60Hz 1Amp on request

Unit 09, 2nd Floor, Prabhadevi Indl Est, Prabhadevi MUMBAI-400025, Tel: - 022-24220297,
E-Mail:- info@campbellelectronics.com / tabm@rediffmail.com, www.campbellelectronics.com