

# Disintegration Tester

## BJ-1/2/3

BJ-1/2/3 Series Disintegration Tester Series are the latest development and production of our company, which is the requisite instrument in detecting disintegration of solid dosage products. The design and technical specs of these testers reached or surpassed the advanced level of this kind of instrument in the world. All our lab instruments are engineered and built in according to USP Specifications. Power source: 220V/50HZ/600W



**BJ-1**



**BJ-2**



**BJ-3**

### Information

- Disintegration series used a MPU to control all operations. The bath liquid can achieve an even temperature.
- Each test station is driven independently. It stops running at the end of a test..
- All components are selected of high quality brands.
- It can be operated conveniently by using keys. Its measurement results are displayed by LED.
- High quality designed and ruggedly built to meet the USP Codes.

SPECIFICATIONS	BJ-1	BJ-2	BJ-3
Nacelles	1	2	3
Nacelle glass tubes	6	12	18
Vol. of Burning Cup	1000 ml		
Nacelle back & forth frequency	30 ~ 32 /min		
Nacelle back & forth range	55 ±1 mm		
Temp. Range	Ambient to -45℃		
Automatic Temp. Accuracy	±0.5 °C		
Power of the Heater	500W		
Continuous Work	> 24 hours		
Power	220V/50HZ/600W		
Dimension (L x W x H)mm	260 * 300 * 420	400 * 300 * 420	540 * 300 * 420
Weight	11kg	16kg	22kg