

Tap density Model AS-100



Tapped Density Analyzer or Tapped Density Meter . AS-100 Tapped Density Tester is a density measurement equipment used to measure the tapped density of powders, as well as granulated or flaked materials by standardized and repeatable procedures

Technical Specifications

- Sample weight: 500.00g
- Tapped sample volume: 250.00ml
- Number of tapping: 0-99999 (set 0 to output bulk density)
- Motor torch: 0.86 N.m.
- Frequency: Maximum 300rpm (adjustable)
- Magnitude: Maximum 14mm (adjustable in integer from 1-14mm)
- Input power: AC 220V (50Hz/60Hz)