



Hotplate Magnetic Stirrer รุ่น SH-3

คุณสมบัติ

- พลังงานความร้อนและความเร็วในการกวนสามารถปรับได้แบบ stepless
- มาพร้อมชุดจับเทอร์โมมิเตอร์แบบแท่งสแตนเลส
- ขนาดแผ่น 170 × 170 มม.
- กำลังมอเตอร์ 20 W
- ทำอุณหภูมิได้สูงสุด 380 องศา
- ความเร็วรอบ 100 ~ 2000 รอบต่อนาที
- รองรับปริมาณตัวอย่าง 5,000 มล.
- รับน้ำหนักตัวอย่างได้ 10 กิโลกรัม
- โหมดการแสดงผลความเร็วแบบสเกล
- โหมดการตั้งค่าแบบปรับหมุน

รายละเอียด

Model	SH-3	
Function	Max. Volume (ml)	5000
	Speed (r p m)	100~2000
	Working time	Continuous
	Minimum revolution (r p m)	100
	Max Temp (°C)	380

Structure	Exterior	Cast aluminum, Surface spraying
	Plate size(mm)	170×170
	Driving mode	Motor
	Motor power (W)	20
	Heating Power(kW)	0.5
	Stir bar	PTFE
Controller	Continuous operation	Yes
	Setting mode	Scale
	Speed display mode	Scale
	Timer	No
	Operation function	No
	Temp. Sensor	No
Specification	Exterior size (W*L*H)(mm)	265×185×190
	Packing size (W*L*H)(mm)	280×210×150
	Plate size(mm)	170×170
	Plate load	10KG
	Power supply(50/60Hz)Current rating	110V or 220V
	NW/GW (kg)	2.2/3.4