



SF6

## Features

- ▶ For flocculation test of water, possible to stir 4 or 6 samples at the same time.
- ▶ The fixed illuminator provides glare-free lighting to provide diffused lighting of samples
- ▶ Life time keeps long due to the excellent & durable brushless DC motor.
- ▶ Speed control maintains 4 or 6 peddle's rotation speed concurrently by feedback control
- ▶ Digital speed indication & built-in timer(99hr 59min)
- ▶ Indicate separately set-in rpm an practical rpm on the bright LED screen and confirm at the same time.
- ▶ SP4 : portable type



SF4

## Specifications

Model		SF6	SF4
No. of places		6	4
Stirring	Capacity	up to 2000ml	
	Motor	Brushless DC motor 60W	
	Speed	20~300rpm	
Speed Controller		Feedback control	
Speed Display		Digital	
Timer		99hr 59min	
Dimensions (W×D×H)		900×250×450mm	320×320×390mm
Electrical Supply		220V 50/60Hz	
Weight		21.5Kg	9Kg
Price			