



## Blaine Air Permeability Apparatus model SBT-127

### I. Purpose

For the air permeability Measuring of cement specimen. The testing conforming with the requirements of GB8074 ,ASTM C204 ,BS4359-2

### II. Technical parameters

Model	SBT-127 (Digital)
Inside diameter of permeability cell	Φ12.7 mm.
Height of sample in the cell	15 mm.
Perforated disk	35
Diameter of the hole	Φ1.0 mm.
Thickness of disk	1.0 mm.
Net weight	≈ 6 kg.
Power	220V/50Hz/25W