

# JAR TESTER

## JT-X-PLC series

JT-4-PLC



JT-6-PLC



Automatic Up&Down



Speed Adjustable



LED Light



PLC Touch Screen

Model	JT-4-PLC	JT-6-PLC
No. of places	4 Position	6 Position
Capacity	500 - 1,000 ml	
Speed	10 - 320 rpm	
Display	Touch Screen	
Timer range	999 min 59 sec	
Rotation	Clockwise and Anticlockwise	
Step	4 Step Setting Timer and Speed	

\*Note : Custom-size is Available.

