



กล่องถุงมือ Acrylic Glove Box YKG Series

Technical Features

Transit box vacuum degree	0~ -0.1Mpa (2.5 Vacuum gauge), pressure maintain period for 12Hours above
Operation box vacuum degree	0~ -0.1Mpa (2.5 Vacuum gauge), pressure maintain period for 12Hours above

Product Parameters

Model	Operation Box (mm)	Thicknes s of Acrylic (mm)	Transit Box (mm)	Vacuum Function
YKG-012B	700*450*500	10mm	240*240*240	It owns vacuum function for operation box, and no vacuum function for transit box
YKG-011B	800*550*600	10mm	240*240*240	
YKG-010B	1000*500*500	10mm	240*240*240	
YKG-009B	1000*700*500	10mm	240*240*240	
YKG-008B	900*600*700	10mm	240*240*240	
YKG-007B	1000*600*700	10mm	240*240*240	
YKG-006B	1000*750*700	10mm	240*240*240	
YKG-005B	1200*600*600	10mm	240*240*240	
YKG-004B	1200*600*700	10mm	240*240*240	
YKG-001B	1094*640*645	10mm	240*240*240	

