

AL Series Halogen Rapid Moisture Analyzer



Features:

Simple operation, no need to install, small size, convenient and practical, wide application range, high measurement accuracy, short detection time.

Technical Specification:

Model	AL-1201E	AL-1201
Accuracy	0.001g	
Max.Capacity	120g	
Readability	0.01%	
Weighing sensor	China	HBM
Temperature range	40-160°C	40-230°C
Heating source	Halogen Lamp	Halogen Lamp 450W
Heating way	Gentle, standard, fast	
Halt mode	Automatic, manual, timed	
Result display	Test temperature, test time, moisture content, solid content, dry weight, dynamic curve, moisture regain, solid content, weight loss, dry weight, test method, sample name	
Tray material	Stainless	
Tray diameter	Φ100mm	Φ90mm
Interface	Standard RS232 communication interface, which can be connected to printers, PCs and other peripheral devices, and USB can be customized	
User management	-	Password protection, multiple users, auto-lock
Display Screen	4.3 inch	5 inch HD
Test Data Storage	100 groups	
Product size	290x200x180mm	370x215x195mm
Weight	NW:2.26kg,GW:3.35kg	NW:3.8kg,GW:4.8kg
Standard configuration	One host machine, one stainless steel sample plate, one tripod stand, one heat collecting plate, one calibration weight, one certificate of conformity, and one instruction manual	

