

Dry Ice Logger

Extra-low Temperature PDF data logger

Measurement range : -80°C to 70 °C (-112°F to 158°F)

- 17,280 samples memory capacity
- Replaceable battery
- Activation button to start logging
- USB interface for PC & MAC communication
- Built-in LED for logging and alarm indication
- User defined Alarm Levels, Sample Rate, and Time Zone settings
- Water and dust resistant IP67



Solutions



Hospital
Environments



Universities &
Research Centers



Pharmaceutical
Manufacturing sites



Pharmacies



Laboratories

www.scilution.co.th | Youtube : SCILUTION

Tel : 02 964 8018 | Line ID : @scilution